

# **ARPENspin Series DNA/RNA Extraction Kits**

Product Catalog

# ARPENspin DNA Extraction Kit Series

## Find the Right DNA Extraction Kit for Your Sample

### ARPENspin Omni Genomic DNA Extraction Kit

Fast genomic DNA isolation from diverse samples in 20 minutes.

**Sample types:** Blood, bacteria, cells, tissue, yeast, swabs

**DNA conc.:** 20–200 ng/μL

**Cat. No.:** ARP-GDES50

**Package:** 50 prep



### ARPENspin Plant Genomic DNA Extraction Kit

Fast genomic DNA isolation from plant tissues, including polysaccharide and polyphenol-rich samples.

**Sample types:** Leaves, roots, stems, soybeans, corn

**DNA conc.:** 20–100 ng/μL

**Cat. No.:** ARP-PGDES50

**Package:** 50 prep



### ARPENspin FFPE Tissue DNA Extraction Kit


High-quality DNA extraction from FFPE and formalin-fixed tissue samples

**Sample types:** FFPE or formalin-fixed tissues

**DNA conc.:** 300–500 ng/μL

**Cat. No.:** ARP-FTDES50

**Package:** 50 prep



# ARPENspin RNA Extraction Kit Series

## Find the Right RNA Extraction Kit for Your Sample

### ARPENspin Total RNA Extraction Kit


High-quality total RNA extraction from diverse sample types in under one hour.

**Sample types:** Blood, tissues, cells, bacteria, fungi

**RNA conc.:** 200–400 ng/μL

**Cat. No.:** ARP-TRES50

**Package:** 50 prep



### ARPENspin Plant Total RNA Extraction Kit

RNA extraction from plant tissues, including polysaccharide and polyphenol-rich samples.

**Sample types:** Leaves, roots, stems, soybeans, corn

**RNA conc.:** 20–100 ng/μL

**Cat. No.:** ARP-PTRES50

**Package:** 50 prep



### ARPENspin FFPE Tissue RNA Extraction Kit


High-quality RNA extraction from FFPE and formalin-fixed tissue samples

**Sample types:** FFPE or formalin-fixed tissues

**RNA conc.:** 20–100 ng/μL

**Cat. No.:** ARP-FTRES50

**Package:** 50 prep



### ARPENspin Fast Soil/Stool DNA/RNA Extraction Kit

Simultaneous DNA/RNA extraction from challenging soil and stool samples with optimized inhibitor removal.

**Sample types:** Soil, stool

**Purity:** OD260/280: 1.7–2.2

**Cat. No.:** ARP-FSSDR50

**Package:** 50 prep



# ARPENspin Series | Built for Reliable Extraction

## Confidence in every sample

ARPENspin extraction kits simplify nucleic acid purification from routine and challenging sample types, supporting clean DNA and RNA for downstream molecular workflows.



## What ARPENspin Delivers



### Safe & Non-Toxic

Reduced hazardous reagent use for safer workflows



### Fast & Simple

Simple spin-column workflows for efficient extraction



### High Purity Output

Pure nucleic acids for PCR, qPCR, RT-PCR, and sequencing.



### Reliable with Complex Samples

Optimized workflows for plant, FFPE, soil, and stool samples

**Fast, easy-to-use, and cost-effective extraction workflows designed to deliver high-quality DNA and RNA with reliable yields.**



## Explore the full ARPENspin extraction kit range

Find detailed product information, specifications, and ordering details for ARPENspin DNA, RNA, and DNA/RNA extraction kits.

Visit [www.arpengen.com](http://www.arpengen.com) to learn more.

[Visit Website](http://www.arpengen.com)